



2024

KOHEZIVNI POVOJI



VELEMEDICA



KOHEZIVNI POVOJ

Raztegljivost: 86,9%, Oprijemljivost: 54 cN/cm, Sila odvijanja: 84 cN/cm

Samosprijemljivi kohezivni elastični bandažni povoj z odlično raztegljivostjo, dobrim sprijemanjem in enostavnim trganjem.

Povoj se ne sprijema s kožo in lasmi. **Brez lateksa.**

Kataloška številka	Dimenzija	Barva	Osnovno pakiranje (SC)	Transportno pakiranje (BOX)
VCB7545-B	7,5cm x 4,5m	Modra	12 kos	16 SC
VCB7545-G	7,5cm x 4,5m	Zelena	12 kos	16 SC
VCB7545-R	7,5cm x 4,5m	Rdeča	12 kos	16 SC
VCB7545-Y	7,5cm x 4,5m	Rumena	12 kos	16 SC

Osnovna šifra: FBZN-01

**možna dobava različnih variacij dimenzij in barv.*